

## Micro Cannula Needle 16G 100mm For Dermal Filler

Our Product Introduction

for more products please visit us on hadermal-filler.com

### Basic Information

- Place of Origin: China
- Brand Name: Blunt Micro Cannula
- Certification: CE
- Model Number: 14g-30g
- Minimum Order Quantity: 30pcs
- Price: 0.8usd/pcs
- Packaging Details: bag
- Delivery Time: 4-7work days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 13333box per week



### Product Specification

- Product Name: Micro Cannula
- Material: Stainless Steel Blade
- Function: Micro Cannula For Dermal Filler Injection
- Application: Medical Disposable Blunt Tip Needle
- Usage: Single-use
- Package: 1pc/Pack,50pcs/box
- Type: Medical Disposable Blunt Tip Needle Micro Cannula
- Keyword: Blunt Tip Needle Micro Cannula
- Highlight: **16G Micro Cannula Needle, Dermal Filler Micro Cannula Needle, 100mm Micro Cannula Needle**



### More Images



### Product Description

#### Medical Disposable Blunt Tip Needle Micro Cannula For Dermal Filler Injection

#### Micro cannula blunt tip needle for beauty injection

#### Specifications:

- 1) Superior quality and complies with international cosmetology and medical requirements.
- 2) Needle made of medical grade stainless steel.
- 3) sterile package

Our Product

- 4) Disposable, for single use only.
- 5) Different size available, special design is available

**What's The Procedure?**

At the source of the filler injection using Microcannula, your Doctor creates a minute entry point in your skin using a precise, sharp needle. This is a vital step since the Microcannula has a blunt tip and is incapable of penetrating the skin by itself. The insertion location is minimal and does not cause substantial pain or bleeding. Then the Doctor skillfully inserts the Microcannula into the small opening and glides it under your skin to reach the sections of filler injection. The Microcannula then slides within the spaces connecting the tissue. It can be moved under the skin without any injuries to tissue. Patients typically do not feel much pain, if any, as the Doctor guides it under their skin.

Place of Origin	China
Disinfecting Type	EOS
Material	Stainless Steel
Safety standard	ISO 13485:2003
Application	Face, Neck, Chest and Body More



Item	Size	Needle Length	Product
Blunt Tip Needle (Micro Cannula)	18G	50mm	
		70mm	
		80mm	
	20G	100mm	
	21G	38mm	
	22G	70mm	
	23G	38mm	
		50mm	
		38mm	
	25G	40mm	
		50mm	
		70mm	
	26G	40mm	
		50mm	
		70mm	
27G	40mm		
	50mm		
30G	50mm		
	40mm		



